



## ISO/134.2 KHZ MICROCHIPS

Contact Us: 800-252-7894 [microchip@akcreunite.org](mailto:microchip@akcreunite.org)

Visit Us Online: [www.akcreunite.org](http://www.akcreunite.org)




### INDIGO ISO MICROCHIPS

- 12 gauge needle
- Microchip dimensions: 11.5mm x 2.12mm
- Ideal for all pets
- Needles are easily separated for quick sharps disposal

### INDI ISO MINICHIPS

- 15 gauge needle
- Microchip dimensions: 8mm x 1.4mm
- Ideal for smaller pets, thin skinned animals with smooth coats, and exotic pets
- Double-click implantation tactile feedback for reliable implantation



### AKC Reunite microchips provide permanent pet identification that lasts a lifetime!

- Unique ID Number beginning with 956
- ISO/134 kHz/International Microchips — detected by ISO/Universal Scanners
- Approved for international and domestic travel
- Includes collar tag with microchip number (optional) and 6 barcode labels
- Sold as singles, 5-count boxes and 25-count boxes - ask about large-volume discounts
- Option to upgrade microchips with prepaid enrollments

### Quality-first, superior technology to protect the pets that enhance and enrich our lives!

- Human-grade materials and stringent manufacturing process
- Bio-compatible glass microchip encapsulation with anti-migration Parylene coating
- Pre-sterilized and pre-loaded microchip implanter
- Extra sharp, silicone-coated, double-bevel needle
- RoHS compliant with no lead or other hazardous substances
- Made by Trovan in Europe

### CUSTOMIZE YOUR PROGRAM

- Choose your preferences in The HUB, your online portal for easy records management.
- Lock as an Alternate Contact to help when a pet is found
- Lock as a Primary Contact to be notified if a pet is found and prohibit enrollment transfer

### EXTRA SAVINGS & PERKS





- Free transfers through THE HUB or your pet management software
- No annual fees or renewal fees
- Dedicated Account Manager
- Product replacement assistance

Always check the chip at [aaha.org/petmicrochiplookup](http://aaha.org/petmicrochiplookup)


# MICROCHIP SCANNERS

Contact Us: 800-252-7894 [microchip@akcreunite.org](mailto:microchip@akcreunite.org)

Visit Us Online: [www.akcreunite.org](http://www.akcreunite.org)

	QUICKSCAN 650	PROSCAN 700	WANDSCAN 900 REACH (Pole)	WANDSCAN 900 PRO (Halo)
<b>UNIVERSAL SCANNERS</b> <i>read all brands of pet microchips</i>				
<b>DIMENSIONS</b>	Width 69.8mm (2.7 in.) Length 125.4mm (4.9 in.) Height 24mm (0.9 in.)	Width 120mm (4.72 in.) Length 238mm (9.37 in.) Height 35mm (1.38 in.)	<b>Reader:</b> See QuickScan 650 <b>Pole:</b> grip 35x100 + 500x20mm +2.3cm diameter	<b>Reader:</b> See QuickScan 650 <b>Pole:</b> grip 35x100 + 400x20 + 125mm halo
<b>WEIGHT</b> <i>with batteries</i>	120gm (4.2 oz)	300gm (10.6 oz)	Pole: 70gm	Pole+ Halo 130gm
<b>BATTERY POWER</b> <i>continuous read</i>	9 V battery Approx. 1 hour	4x 1.5V AA battery alkaline/ NiMH rechargeable Approx. 3 hours	9 V battery Approx. 1 hour	9 V battery Approx. 1 hour
<b>INTERFACE</b>	n/a	mini n-USB port	n/a	n/a
<b>TYPICAL READING TIME</b>	Trovan™: 24ms ISO FDX-B : 100ms Destron: 100ms Fecava: 100ms	Trovan™: 24ms ISO FDX-B : 100ms Destron: 100ms Fecava: 100ms	Trovan™: 24ms ISO FDX-B: 100ms Destron: 100ms Fecava: 100ms	Trovan™: 24ms ISO FDX-B: 100ms Destron: 100ms Fecava: 100ms
<b>READING DISTANCE</b>	FDX-B up to 50mm (1.97in.) ID-100 up to 50mm (1.97 in.) Destron clones up to 25mm (1 in.) AVID encrypted up to 30 mm (1.18 in.)	FDX-B up to 80mm (3.2 in.) ID-100 up to 80mm (3.2 in.) Destron clones up to 70mm (2.7 in.) AVID encrypted up to 60 mm (2.4 in.)	FDX-B up to 50mm (1.97in.) ID-100 up to 45mm (1.97 in.) Destron clones up to 25mm (1 in.) AVID encrypted up to 30 mm (1.18 in.)	FDX-B up to 65mm (2.6in.) ID-100 up to 60mm (2.5 in.) Destron clones up to 60mm (2.7 in.) AVID encrypted up to 50 mm (2.4 in.)

## ISO-SCAN 134 - Only reads 134.2 khz/ISO microchips

	<b>DIMENSIONS</b>	Width 70mm (2.75 in.)   Length 154mm (6.06 in.) Height 18mm (0.7 in.)
	<b>WEIGHT</b>	174gm (6.1 oz) - <i>with batteries</i>
	<b>BATTERY POWER</b>	9 V battery   Approx 10 hours ( <i>continuous read</i> )
	<b>INTERFACE</b>	Bluetooth
	<b>READING PERFORMANCE</b>	Typical read time: 400ms   Reading distance: up to 90mm (3.5 in.)

Based on AKC Reunite recovery statistics, more than 90% of microchipped pets in the US have 134.2 khz/ISO microchips.